

PRODUCT

Instrument and accessories



BRAND 18 years of IVD experience; intellectual property rights; the only automatic POCT of its kind.

Technology and Innovation:
Through international collaboration, we have optimized this control system for accurate performance.
















Patents:
52 international patents and 12 invention patents.


Social responsibility :
Sponsor of Taihang Mountain water source infrastructure; participation in relief of Sichuan Earthquake-stricken area; Silver medal of Chinese Youth Volunteer service.

Our Core value and Facilities:
We do our utmost to achieve: simplicity, integrity, diligence, and outperformance.
Our production facility has grown from 500m² to a 10000m² plant.



TEST LIST

-  HCY
-  D-Dimer
-  Anti-CCP
-  PCT
-  NGAL
-  H-FABP
-  cTnI
-  HsCRP
-  HbA1c
-  mAlb
-  Myo
-  ASO
-  RF
-  Cys-C
-  ADA



R&D: 3938 Trust Way, Hayward, CA 94545, USA
Production base: High-Tech Development Zone, Hebei, China
Service Hotline: 400-0191-606

Telephone: +86 311 83832628 | Fax: +86 311 83859087
Email: 4000191606@hipro.us | Website: <http://www.hipro.us>

A3 [INNOVATIVE DESIGN]
Total Laboratory Automation
ECO platform




Hipro[®]
Hipro Biotechnologies Co., Ltd.



400-0191-606
www.hipro.us www.hiprochina.com

OPERATION

Accuracy and simplicity



3 Key Technologies INNOVATION

Automated platform

Multiple test methodologies integration
automatic laboratory platform
International reagent

Full-function module

Single module capable of all functions
combination of 1 ~ 48 modules
Appropriate for laboratories of all sizes
Scalability for cost savings

3-level Calibration system

Physical calibration by standard liquid traced to ANSI
Remote standard curve calibration
System calibration by custom quality control supplied by Bio-Rad

DETAILS

convenient, fast, and meticulous

Automatic sample adsorption
Quantitative tips remove the need for pipettes.

Preloaded reagents
Simplify procedures with accuracy and ease.

Integrated Sample Collector
Eliminates impacts of anticoagulants.

Temperature Control System
Unaffected by environmental fluctuations.
Wind tunnel design allows quiet and quick response.

Scalable Modular Design
Appropriate for laboratories of all sizes, and outdoor use.

QUALITY

RSP-USA is going through JCTLM review. Hipro is to become a JCTLM laboratory. The Hipro assays and control system are certified by ANSI in Aug, 2013. Accredited by (CNAS) ISO17025 and ISO15195 can be expected soon.

TEST MENU

Multiple available reagents

CATEGORY	TEST	CATEGORY	TEST
Infection	PCT	Rheumatoid	ASO
	Hs-CRP		CRP
	SAA		Anti-CCP
Cardiovascular	H-FABP	Renal Function	RF
	MYO		CYS-C
	cTnI		mAlb
	HCY		NGAL
	NT-ProBNP		RBP
	Lp(a)		β 2-MG
Hormones	Lp-PLA2	Gastric function	α 1-MG
	CK-MB		HP
	HCG		PGI
Coagulation	LH	Autoimmunity	PGII
	D-Dimer		IgA
Diabetes	HbA1C		Liver Function
	ADA	IgE	
Liver Function	TRF	Microelement	
	PA		25-OH-VD
	HPT		More to come

PARAMETERS

Data interfaces: LAN, LIS/HIS, RS232, Bluetooth, Wi-Fi, USB, SD
Power: AC 100~265V
Volume: L xW x H: 385 x 305 x 353mm
Weight: 15Kg